

HEXAR is a leader in the wearable exoskeleton robot technology





Better Life with Robot 2022 09





본사 : (15588) 경기도 안산시 상록구 한양대학로 55 산학연협력관 301~303호 Tel 031 520 1550 Fax 031 520 1560 www.hexarhc.com 강남지점 : (06170) 서울특별시 강남구 영동대로106길 23(삼성동) 중앙빌딩 1층







Integrated active rehabilitation training equipment for recovering of knee joint exercise, muscular strength measurement and muscle reconstruction



Resilion K30

Product Overview

Integrated active rehabilitation training equipment for recovering of knee joint exercise, muscular strength measurement and muscle reconstruction

Main Function


- Isokinetic rehabilitation exercise for patients with sarcopenia and surgery
- Multimodal comprising of Isokinetic, Isotonic, Isometric, CPM functions
- ROM measurement and isokinetic motion test
- Graphic analysis of motion
- Easy left/right switching
- User registration, exercise data history management and print out

Specification

| mm (D) × 1 |
|------------|
| |
| |
| Isotonic , |
| |
| |

Picture





| Π | ٨ |
|---|---|
| U | H |

1275mm (H)

/ Isometric



Specification is subject to change without notice.